

COLD CABINET: REFRIGERATOR

IG-RF



CONTACT US

info@igenels.com
www.igenels.com

ABOUT THE COMPANY

iGene Labserve Pvt. Ltd. is gaining recognition by offering resilient, innovative solutions in laboratory instrumentation across healthcare, genomics, drug discovery, biopharma, and food & beverage sectors. We strive to enhance lab efficacy and reduce challenges through advanced technologies and a diverse product portfolio tailored to our customers' needs.

DESCRIPTION

The IG-RF refrigerator is designed to maintain products at an optimal temperature range of 2–8°C ensuring ideal storage conditions and extended shelf life. It features an manual-defrost system for hassle-free operation, making it a reliable solution for temperature-sensitive materials.

KEY FEATURES

Manual Defrosting:

- Requires periodic user-initiated defrosting to maintain optimal cooling efficiency and prevent ice build-up.

Interior Lighting:

- Includes a built-in interior light for clear visibility and easy access to stored items, even in low-light conditions.

Precise Temperature Control (2–8°C):

- Designed to maintain a stable internal temperature between 2°C to 8°C ideal for storing temperature-sensitive materials and extending product shelf life.

Adjustable Shelving:

- Available with 5 to 10 shelves, depending on the model capacity, offering flexible storage options to suit varying laboratory or clinical needs.

Durable Construction:

- Made from high-quality die-cast plastic, ensuring durability, lightweight handling, and long-lasting performance in demanding laboratory environments.

Customizable Options:

- Customization is available to meet specific user requirements, ensuring the refrigerator aligns perfectly with operational needs.

TECHNICAL SPECIFICATION

Model	IG-RF
Gross volume available	200L to 600L
Temperature range	2°C - 8°C
Power	225 W
Inner chamber	Die-cast plastic body
External cabinet	Die-cast plastic body
No. of door	1 Door
No. of shelves	5 Shelves
Interior light	Yes
Wheels/castors	Yes
Refrigerant	R134a
No. of compressor	1 Compressor
Optional	Customization is available



CONNECT WITH US.

ADDRESS

16/2, BLOCK 16, ASHOK NAGAR,
NEW DELHI-110018

TELEPHONE

1800-572-0603

 **EMAIL**
info@igenels.com

 **WEBSITE**
<https://www.igenels.com>



GeM
Government
e Marketplace



एन एस आई सी
NSIC
ISO 9001 : 2008

