

Refrigerator

IG-RF Series



IGENE LABSERVE
Innovative • Interactive • Intuitive

FEATURES & SPECIFICATIONS

DESCRIPTION

IG-RF series is suitable for keeping the products at optimum temperature. The refrigerators are equipped with auto defrosting device for ease of operation. This can maintain a temperature ranging from 2-8°C and hence can provide optimum conditions to store the material and enhance the shelf life of the products stored in it.

FEATURES

- ❖ Automatic defrosting
- ❖ Comes with an interior light
- ❖ Can maintain a temperature ranging from 2-8°C
- ❖ Comes with 5-10 shelves depending on the capacity of the refrigerator

SPECIFICATIONS

Model number	IG-350RF
Gross Volume (Litres)	350
Temperature Range	2-8°C
Power (W)	225
Lock	Yes
Defrost	Automatic
No. of Lids/Doors	1
No. of shelves	5
Interior light	Yes
Wheels/castors	Yes
No. of compressors	1
Refrigerant	R134a

Optional Features:

Non-transparent doors can be provided as per the requirement of the customer.