


Horizontal Electrophoresis

IG-HE Series



IGENE LABSERVE

Innovative • Interactive • Intuitive

FEATURES & SPECIFICATIONS

FEATURES

- ❖ Made of high-quality Polycarbonate, exquisite and durable.
- ❖ Transparent, convenient for observation.
- ❖ Easy and simple for the end user to use.
- ❖ Special gel device for easy and fast gel casting.
- ❖ High resolution and simple separation.
- ❖ Anti-high temperature, in deformable.
- ❖ Withdrawable electrodes convenient for maintenance.
- ❖ Color coded electrodes for easy operation.
- ❖ One gel casting device can cast different size of gels.



SPECIFICATIONS

Model	IG-20HE	IG-25HE
Gel dimension	20x20 cm & 20x25 cm	25x25 cm & 33x25 cm
Sample combs	4(1.5, 1.5) & 4 (1mm)	4(1.5) & 4 (1mm)
Comb thickness	1.5 mm, 1mm	1.5mm & 1mm
Buffer volume	2.5 L	4L
Max. no. sample	60	104
Size (LxWxH)	43x23.5x7.5	53x28x8.5

