


# Multitube Vortex Mixer IG2-500



**IGENE LABSERVE**

Innovative • Interactive • Intuitive

# FEATURES & SPECIFICATION

## FEATURES

- ❖ Maximum process of 50 samples at the same time to make the experiment more convenient and efficient.
- ❖ Microprocessor control, simplistic designed appearance, LED display both speed and time.
- ❖ Built-in short mix and timing operation modes. Stable and silent operation.
- ❖ Large orbital vortex, effective performance. Various types of foam tube frames and pads for option.



# SPECIFICATIONS

<b>Model</b>	<b>IG2-500</b>
<b>Speed Range</b>	500 ~2500rpm
<b>Speed Accuracy</b>	±10rpm
<b>Shaking Orbit</b>	4mm
<b>Time Range</b>	1s ~9999min, short mix, impulse pattern, continuous
<b>Interval timing between impulse</b>	1 ~10s
<b>Timing set for impulse working</b>	1s ~99min59S
<b>Max Load</b>	5kgs
<b>Cover Plate Size</b>	184 x 311mm
<b>Voltage</b>	AC 100 ~230V, 50/60Hz
<b>Power</b>	60W
<b>Fuse</b>	250V, 1A, Φ5 x 20
<b>Dimension (WxDxH)</b>	250 x 426 x 480mm
<b>Net Weight</b>	15kgs

# ACCESSORIES

<b>Cat #</b>	<b>Description</b>	<b>Dimension</b>
<b>DMT-A</b>	foam tube frame Φ10mm x 50	245 x 132 x 45mm
<b>DMT-B</b>	foam tube frame Φ12mm x 50	245 x 132 x 45mm
<b>DMT-C</b>	foam tube frame Φ13mm x 50	245 x 132 x 45mm
<b>DMT-D</b>	foam tube frame Φ16mm x 50 (suit for 15mltube)	245 x 132 x 45mm
<b>DMT-E</b>	foam tube frame Φ25mm x 15	245 x 132 x 45mm
<b>DMT-F</b>	foam tube frame Φ29mm x 15 (suit for 50mltube)	245 x 132 x 45mm
<b>DMT-G</b>	foam tube frame Φ38mm x 10 (suit for 100mltube)	245 x 132 x 45mm

