

Vortex Mixer

IG-VMSSMini



IGENE LABSERVE

Innovative • Interactive • Intuitive

FEATURES & SPECIFICATIONS

DESCRIPTION

IG-VMSMini vortex mixer is a touch-to-mix (short mix) device that works with a range of tubes and centrifugal tubes. Its tiny design makes it appealing and welcoming, and its touch-to-mix function makes it simple to use and practical. With its prominent mixing effect, it becomes a vital experimental instrument.

FEATURES

- ❖ Maximum diameter of the tube is upto 30mm, suits for mixing tubes/centrifugal tubes.
- ❖ Outer shell and touch concave are made of plastic which is resistant to chemical attack.
- ❖ Thick base, makes operation steady and reliable.
- ❖ 12V low voltage power adapter built in, smart and portable.

SPECIFICATIONS

Model number	IG-VMSMini
Speed Range	0-4000rpm
Switch Mode	Jog Switch + Knob adjust
Shaking orbit	4.5mm (circle)
Max sample capacity	50ml
Power Input / Output	1.2W / 0.8W
Voltage	DC12V 1A
Power	12W
Dimensions	W.110 x D.110 x H.86mm
Net Weight	1.1 Kgs

