

# Mini Centrifuge

**IG-1010**



**IGENE LABSERVE**

Innovative • Interactive • Intuitive

# FEATURES & SPECIFICATIONS

## DESCRIPTION & FEATURES

- ❖ Large radius rotor with 8 tube holes, centrifugal force is 1.5 times bigger than ordinary 6-hc rotor at the same rolling head.
- ❖ The cover press-button is outward designed with start/stop function. Press the button, it will open automatically to 95°, which is very convenient for one handed operation.
- ❖ The rotor is clip-on fixed. Very convenient for rotors exchanging yet fixed solid reliable.
- ❖ LED display speed and timing.
- ❖ Working is silent and stable. Maintenance free motor, long life and safe.
- ❖ Smart appearance. Multi-purpose to meet different experimental demands.

## SPECIFICATIONS

<b>Max. Speed</b>	<b>3000-10000rpm</b>
<b>Max. RCF</b>	700×g
<b>Rotor Capacity</b>	8x2.0ml / 0.5ml centrifuge tubes (8x0.5ml & 8x0.2ml) 2 PCR strips or 2x8x0.2ml tubes
<b>Driving Motor</b>	DC motor
<b>Power</b>	AC110-220V/50Hz/60Hz
<b>Time Range</b>	1s-999min or continuous
<b>Noise level</b>	≤45dB
<b>Dimensions [D × W × H]</b>	176 × 156 × 121mm
<b>Weight</b>	1.5kg

