

Vacuum Pump

IG-VP610



IGENE LABSERVE

Innovative • Interactive • Intuitive

FEATURES & SPECIFICATIONS

FEATURES

- ✦ Superior chemical resistance**
All wetted parts of pumps are made of PTFE and FFKM material, making them resistant to most of corrosive gases.
- ✦ Start against vacuum**
designed to start against vacuum, are ideal for various kinds of laboratory instruments which need to stop and restart pump repeatedly during operating process.
- ✦ Oil-free & care-free**
Vacuum pumps are driven by diaphragm that needs no lubricant, making it maintenance-free and oil-mist-free.
- ✦ Quieter & more stable**
Motor-driven transmission and vibration-proof assembly, making pumps much quieter and more stable.
- ✦ Overheat protection**
All pumps are equipped with overheat protection device that shuts down when overheated and resumes once cooled down, making the pumps safer and long-lasting.

SPECIFICATIONS

Model	IG- VP610
Max power	220 W
Max current	2A
Max vacuum	7mbar abs.
Max flow rate	34 L/min
Noise	60 dB
Hose barb	ID 10
Net weight	13.3 kg
Dimension (LxWxH)	35x16x23 cm